



Digital X-ray top set, 4000-/5000-SERIES tables

BF101



Removable, radiolucent X-ray tops mount onto surgical tables. For use with Hook-and-loop table pads. For use with table manufactured AFTER November 30, 2015.

Head 4.8 lbs (2kg) / Back 10 lbs (4,54kg)

Seat 7.6 lbs (3,45kg) / Leg 10.2 lbs (5kg).

Dimensions:

Head Section 20" (510mm) (L) x 9.65" (240mm) (W)

Back Section 21.06" (540mm) (L) x 9.65" (240mm) (W)

Seat Section 17.72" (45mm) (L) x 9.65" (240mm) (W)

Leg Section 23.43" (590mm) (L) x 9.65" (240mm) (W)

COMPATIBLE WITH

4085 - 5085 SRT - 5085 HARMONY - 4095 - 5095

